


# Digital Camera Viscam



## Product Features

- No additional capture card needed
- Directly connect to PC's USB 2.0 port
- Can be connect directly to notebook or PC with USB 2.0
- Requires no external power supply to camera
- View live high resolution, uncompressed video on PC
- Capture high resolution image
- Support Microsoft Direct Show - WDM, Twain, and VFW drivers
- C-mount, Tripod mounting socket
- Measuring & Image Processing S/W also include
- Support Window XP / Window 7 / Vista

Model	Viscam1320	Viscam3010
Type	1/2" CMOS	1/2" CMOS
Pixels	1.3 Mega pixels	3.2 Mega pixels
Resolution	1280 x 1024	2048 x 1536
Frame rate	15	6