



บริษัท ออปติค้า จำกัด **OPTICA Co.,Ltd.**

9/231-232 หมู่ 5 ถ. ศรีนครินทร์ ต. ลำโพงเหนือ อ. เมือง จ. สมุทรปราการ 10270 โทร. (662) 385-7357-8 โทรสาร (662) 385-7106
9/231-232 Moo 5 Srinakarin RD., T.Samrongnue, A. Muang, Samutprakarn 10270 Tel. (662) 385-7357-8 Fax. (662) 385-7106

E-mail : optica@opticascope.com

Optech Fixed Stereo Microscope Model SFM

Product Features

- Viewing head : 45° inclined, 360° rotates, Binocular head or Trinocular head. The inter-papillary distance adjustable between 54-75mm; both diopter adjustment of range ± 5 mm.
- Eyepiece : Wide field, high-eyepoint eyepiece, It can observe image of observation specimen with glasses.
- Objective : fixed objective magnification 1x, 2x, 3x, 4x
- Working distance : 100mm.
- Illumination : Built-in and Built out illumination system.
- Accessories : There are many accessories to choose.



SFM-T



SFM-B

Model	SFM1020	SFM1030	SFM2040
Eyepiece	WF10x/20	WF10x/20	WF10x/20
Objective	1x ,2x	1x ,3x	2x ,4x
Magnification	10x ,20x	10x ,30x	20x ,40x
Working distance (mm.)	100	100	100

Website:<http://www.opticascope.com>