



## Illuminator ALPHA NEXT

### Product Features

- Streamline facade
- New design for easy replacing lamp
- Built-in Blue filter
- Brightness > 23000LU



**ALPHA Next**

The newest addition to our line of illumination equipment is the Alpha NEXT. It features “cold” illumination through filtration. The cold light source has high intensity fiber optic illumination and accessories that allow you to position your light for best possible viewing. The 150 watt lamp is equipped with a variable control. The ALPHA NEXT provide the best illumination for high brightness & low temperature with light weight.

<b>POWER INPUT</b>		220V
<b>POWER OUPUT</b>		5V-21V
<b>POWER CONSUMPTION</b>		150W
<b>DIMENSIONS</b>		270mm×150mm×130mm
<b>FILTER</b>		Blue
<b>REFLECTOR LAMP</b>		21V/150W
<b>FIBER ADAPTER</b>		F20
<b>OPTION</b>	<b>Twin-Arm Fibre</b>	Ring Fibre
<b>SPECIFICATION</b>	<b>2- F5X500mm</b>	F62XF8X550mm
<b>OPTION: Microscope Adapter for Ring Fiber</b>		
<b>ADAPTER SPECIFICATION</b>		M44X0.75
		M48X0.75
		M50X0.75
		M54X0.75